

NIRS Analyzer PRO – FiberSystem

A629281120


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery A629281120

Qt.	Order no.	Description
1 PCS	6.7440.480	NIRS Process, I/O Communication box
		I/O Communication box for process instruments with 24 analog and 24 digital channels

A circular icon containing a grey rectangular box with a white circle in the center, representing the I/O Communication box. The icon is crossed out with a diagonal line.

Optional accessories

Order no.	Description
A660044277	NIRS Analyzer PRO Lamp 12V/20W
	Replacement lamp for the NIRS Analyzer Pro (12V/20W)

